





This guide may be used in deciding which wound dressing to use in conjunction with V.A.C. VERAFLOR<sup>TM</sup> Therapy

	V.A.C. VERAFLOR <sup>TM</sup> Dressings	V.A.C. VERAFLOR <sup>TM</sup> Large Dressing	V.A.C. VERAFLOR CLEANSE <sup>TM</sup> Dressing	V.A.C. VERAFLOR CLEANSE CHOICE <sup>TM</sup> Dressing
<b>Wound Types</b>				
<b>Wound Characteristics</b>	Open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Large open wounds, including wounds with shallow undermining or tunnel areas where the distal aspect is visible	Cavity wounds or wounds with complex geometries	Wounds with thick fibrinous exudate, slough, infectious material and other wound bioburden
<b>Key Goal(s) of Therapy</b>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR<sup>TM</sup> Therapy, to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Generation of robust granulation tissue</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR<sup>TM</sup> Therapy, to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Generation of robust granulation tissue in large wounds</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR<sup>TM</sup> Therapy, to initiate active therapy and to help facilitate the removal of infectious material and other wound bioburden.</li> <li>Easy application into tunneling and undermining</li> </ul>	<ul style="list-style-type: none"> <li>When used in conjunction with V.A.C. VERAFLOR<sup>TM</sup> Therapy, to initiate active therapy and to help facilitate the removal of infectious material such as thick fibrinous exudate, slough, and other wound bioburden.</li> <li>Subject to clinical assessment, this combination may be considered when surgical debridement is not possible or not appropriate</li> </ul>

	V.A.C. VERAFLU™ Dressings	V.A.C. VERAFLU™ Large Dressing	V.A.C. VERAFLU CLEANSE™ Dressing	V.A.C. VERAFLU CLEANSE CHOICE™ Dressing
<b>Shape</b>	Spiral-cut foam	Block foam pre-slit into two layers	Tubular shape	Block foam pre-slit into three layers
<b>Application Characteristics</b>	<p>Easy application:</p> <ul style="list-style-type: none"> <li>• Size without scissors</li> <li>• Precut area for pad application when used for bridging</li> <li>• Single pad application</li> </ul>	<p>Easy application:</p> <ul style="list-style-type: none"> <li>• Ideal for large surface areas with shallow depths</li> <li>• Provided with V.A.C. VERAT.R.A.C. DUO™ Tube set for extended surface area coverage</li> </ul>	<p>Application flexibility:</p> <ul style="list-style-type: none"> <li>• Ideal for addressing wounds with complex geometries (eg, tunnels, undermining and circumferential applications)</li> <li>• Single pad application</li> </ul>	<p>Easy application with flexibility:</p> <ul style="list-style-type: none"> <li>• Thin layers for improved conformability</li> <li>• Multiple layers provide application options for wounds with varying depths</li> <li>• Ideal for dirty wounds needing active therapy</li> <li>• Single pad application</li> </ul>
<b>Dressing Characteristics</b>	<p><b>Small:</b> 7.7 x 11.2 x 1.75cm per piece (2 pieces spiral)</p> <ul style="list-style-type: none"> <li>• Each piece 24cm uncoiled</li> <li>• Each piece 18cm long bridge</li> </ul> <p>Fill Assist Set Point Range: 12-80ml (1 piece); 26-160ml (2 pieces)</p> <p><b>Medium:</b> 14.7 x 17.4 x 1.75cm per piece (2 pieces spiral)</p> <ul style="list-style-type: none"> <li>• Each piece 80cm uncoiled</li> <li>• Each piece 74cm long bridge</li> </ul> <p>Fill Assist Set Point Range: 38-240ml (1 piece); 80-450ml (2 pieces)</p>	<p>25.0 x 15.0 x 1.6 cm per piece (oval - 2 pieces)</p> <p>Fill Assist Set Point Range: 55-350ml (1 piece); 110-500ml (2 pieces)</p>	<p>1 each coiled rod (24" rod; 1.25" diameter)</p> <p>Fill Assist Set Point Range: 65-130ml (1/2 rod); 130-250ml (full rod)</p>	<p>18.0cm x 12.5cm x Varying Thicknesses (0.8cm – 3.2cm)</p> <p><b>3 Layers</b></p> <p><b>1.6cm Cover Layer</b></p> <ul style="list-style-type: none"> <li>• Fill Assist Set Point Range: 85-160ml</li> </ul> <p><b>0.8cm Cover Layer</b></p> <ul style="list-style-type: none"> <li>• Fill Assist Set Point Range: 42-80ml</li> </ul> <p><b>0.8cm Contact Layer</b> (1.0cm circular holes; 0.5cm spacing)</p> <ul style="list-style-type: none"> <li>• Fill Assist Set Point Range: 24-46ml</li> </ul>

**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. This material is intended for healthcare professionals.**